



LaunchWorks

Letting You Focus on Science 

VIRAL TRANSPORT MEDIA (VTM) ROOM TEMPERATURE STORAGE

Our Viral Transport Media is listed per FDA guidance, “Enforcement Policy for Viral Transport Media during the Coronavirus 2019 (COVID-19) Public Health Emergency.” The VTM is manufactured in the LaunchWorks CDMO facility which is registered, cMGP compliant, and ISO 13485:2016 certified.

Call or email us to learn more! 978-338-3045 or info@launchworksCDMO.com

Launchworks VTM Recipe:

- Fetal Bovine Serum (sterile, heatinactivated, USDA approved)
- Hanks Balanced Salt Solution(no phenol red)
- Gentamicin
- Amphotericin B
- Physiological pH range
- Non-inactivating VTM (no guanidine or guanidinium)

Media Specifications

- Microbial Sterility Tested
- DNase and Protease free Tested



LaunchWorks

Letting You Focus on Science 



123 Brimbal Ave.
Beverly, MA 01915
Phone: 978-338 3045
LaunchWorksCDMO.com

Viral Transport Media (VTM)

Product	VTM Filled Tubes	1L VTM Filled Bottles	bulk Package Kit	Standard Kit	Special Kit
Part Number	LWBVTMCT03	LWBVTMCB1L	LWBVTMCBK	LWBVTMC01	LWBVTMCKSP
Description	3 mL of VTM in IATA RATED 12 mL Tube	1L Bottles	100 labeled tubes & swabs, biohazard bags in bulk	Labeled tube in biohazard	Any combination of additional parts
Configuration	Case of 50 or 500	One (1) bottle or case of six (6)	Case of 100	Cases of 100 or greater	Cases of 100 or greater
Options	Patient Label	None	None	None	None

Storage and Transport Temperature: 2 – 30° Celsius, 24 month storage

VTM Specimen Kit

- 1 tube with 3 mLs of sterile LaunchWorks VTM
- Biohazard bag - tube labeled per CDC instructions
- Absorbent pad

Options:

- Nylon-flocked NP swab
- Barcode or QP-coded labels (up to 3)
- Patient Card
- UN 3373 Return Box

Swab Specifications

- Nasopharyngeal swab
- Individually packaged
- Nylon flocked
- 155mm, 80 mm breakpoint
- FDA-registered Manufacturer
- Ethylene Oxide Sterilization
- 2 year shelf-life
- CE-marked

